



# TRILLIUM DIAGNOSTICS, LLC

INNOVATIVE DIAGNOSTICS FOR CLINICAL CYTOMETRY AND LABORATORY HEMATOLOGY

TRILLIUM DIAGNOSTICS

## CD64 (clone 22)


## ASR

Clone	22
Species	Mouse (anti-human)
Source	Cell culture fluid
Specificity	<p>CD64 reacts with the Fc gamma receptor 1 (Fc<math>\gamma</math>R1), a glycoprotein with a molecular weight of 72 kDa. The antigen is a high-affinity receptor for IgG and is able to bind human IgG<sub>1</sub> and IgG<sub>3</sub>. The antigen is expressed on macrophages, monocytes, and interferon <math>\gamma</math> (IFN-<math>\gamma</math>) and Granulocyte Colony Stimulating Factor (G-CSF) stimulated neutrophils.</p> <p>This antibody binds to a CD64 epitope distinct from the ligand-binding site and from the clone 32.2 antibody to CD64 epitope. It shows especially high affinity binding to human mononuclear phagocytes and polymorphonuclear leukocytes exposed to IFN-<math>\gamma</math> or G-CSF. Its binding is not blocked in the presence of human IgG or immune complexes.</p>
Immunoglobulin Class	IgG <sub>1</sub> $\kappa$
Applications	<ul style="list-style-type: none"><li>• ELISA</li><li>• Flow Cytometry</li><li>• Fluorescent Microscopy</li><li>• Western blot, Immunohistochemistry</li></ul>
Availability	<p><b>FITC</b> – Cat. No. CD64-22F 50 <math>\mu</math>g/mL ■ 50 <math>\mu</math>g in 1.0mL of PBS containing 0.2% BSA and 0.09% NaN<sub>3</sub> Working dilution is best determined by user.</p> <p><b>PE</b> – Cat. No. CD64-22P 25 <math>\mu</math>g/mL ■ 25 <math>\mu</math>g in 1.0mL of PBS containing 0.2% BSA and 0.09% NaN<sub>3</sub> Working dilution is best determined by user.</p> <p><b>PE-Cy5</b> – Cat. No. CD64-22C 50 <math>\mu</math>g/mL ■ 50 <math>\mu</math>g in 1.0mL of PBS containing 0.2% BSA and 0.09% NaN<sub>3</sub> Working dilution is best determined by user.</p> <p><b>Purified</b> - Cat. No. CD64-22U 100 <math>\mu</math>g/mL ■ 100 <math>\mu</math>g in 1.0mL of PBS containing 0.2% BSA and 0.09% NaN<sub>3</sub> Working dilution is best determined by user.</p>
Handling & Storage	Store vials upright, tightly capped, at 2-8°C when not in use. Unopened vials are stable until the expiration date indicated on each vial. Avoid unnecessary cycles of warming and cooling. Protect product from temperatures above 30° C and from prolonged time at room temperature (18-26°C). Conjugated antibodies should not be frozen and should be protected from light.
Warning	Sodium azide (NaN <sub>3</sub> ) is a toxic and dangerous compound when combined with acids or metals. Solutions containing sodium azide should be disposed of properly.

*Analyte Specific Reagent ■ Analytical and performance characteristics are not established.*

# CD64 (clone 22)

## Characterization



This antibody has been found to be effective for direct immunofluorescence staining of human myeloid cells for flow cytometric analysis when used at the concentration of 10 $\mu$ l / 10<sup>6</sup> cells. However each investigator should titrate the antibody in their application to determine the optimal per-test amount.

The performance and specificity of this antibody have been tested using Trillium's in-house quality control methods. Manufacturing of the antibody preparation is done using good manufacturing production guidelines.

## Warranty

This product is made in the EU and is warranted to conform to the labeled specifications. There are no warranties expressed or implied that extend beyond the labeled product description. Neither Trillium Diagnostics nor IQ Products will be held liable for any damage to person or property, or any economic loss related to the use of this product. Trillium Diagnostics' or IQ Products' sole liability is limited to replacement or refund of defective product.

## Technical Support

European distribution and technical support provided by IQ Products • +31 (0)50 57 57 000 • orders@iqproducts.nl

## References

- 1 Akerley WL, Guyre PM, Davis BH (1991) Neutrophil activation through high-affinity Fc gamma receptor using a monomeric antibody with unique properties. *Blood* **77**, 607-615.
- 2 Pizza F, Davis BH, Hendrickson S, Mitchell J, Pace J, Bigelow NC, DiLauro P, Nagieri T (1996) Adaptation to eccentric exercise: effect on CD64 and CD11b/CD18 expression. *J Applied Physiology* **80**, 47-55.
- 3 Schiff D, Rae J, Martin T, Davis BH, Curnutte J (1997) Increased phagocyte CD64 expression and improved Fc-receptor mediated phagocytosis following in vivo recombinant human interferon-g treatment of normal human subjects. *Blood* **90**, 2987-94.
- 4 Pan LY, Mendel DB, Zurlo J, Guyre PM (1990) Regulation of the steady state level of Fc gamma RI mRNA by IFN-gamma and dexamethasone in human monocytes, neutrophils, and U-937 cells. *J Immunol* **145**, 267-75.
- 5 Petroni KC, Shen L, Guyre PM (1988) Modulation of human polymorphonuclear leukocyte IgG Fc receptors and Fc receptor-mediated functions by IFN-gamma and glucocorticoids. *J Immunol* **140**, 3467-72.
- 6 Davis, BH, Bigelow, NC, Curnutte JT, Ornvold, K (1995) Neutrophil CD64 expression: Potential diagnostic indicator of acute inflammation and therapeutic monitor of interferon-g therapy. *Lab Hematol* **1**, 3-12.
- 7 Guyre, PM, Campbell, AS, Kniffin, WD, Fanger, MW (1990) Monocytes and polymorphonuclear neutrophils of patients with streptococcal pharyngitis express increased numbers of type I IgG Fc receptors. *J Clin Invest* **86**, 1892-1986.
- 8 Davis BH (2005) Improved diagnostic approaches to infection/sepsis detection. *Expert Rev Mol Diagn* **5**:193 – 207.
- 9 Herra CM, Keane CT, Whelan A (1996) Increased expression of Fc gamma receptors on neutrophils and monocytes may reflect ongoing bacterial infection. *Journal of Medical Microbiology* **44**, 135-40.
- 10 Leino L, Sorvajarvi K, Katajisto J, Laine M, Lilius E, Pelliniemi T, Rajamaki A, Silvoniemi P, Nikoskelainen J (1997) Febrile infection changes the expression of IgG Fc receptors and complement receptors in human neutrophils in vivo. *Clin Exp Immunol* **107**, 37-43.
- 11 Davis BH and Bigelow NC (2005) Comparison of neutrophil CD64 expression, manual myeloid immaturity counts, and automated hematology analyzer flags as indicators of infection or sepsis. *Lab Hematol* **11**:137 – 147.
- 12 US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910. 1030, Occupational Exposure to Bloodborne Pathogens; Final Rule. Federal Register 235:64175-82. 1991
- 13 US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. HHS Publication (NIH) 93-8395. Washington: US Government Printing Office. 1993

*Analyte Specific Reagent* ■ *Analytical and performance characteristics are not established.*