

# Gentech

**DIAGNOSTICS**



 **Gentech Diagnostics P Ltd**, sole distributor for  **IQProducts**  
bright fluorescence

 **IQProducts**  
bright fluorescence

**Price Catalogue - FY 2010**

234-A, First Floor, Chandi Plaza Sant Nagar, East of Kailash, New Delhi-110065, India  
Tel : +91-11-26445263, 26429843, 46517353; Fax : +91-11-26445269  
E-mail : [gentech@gentechin.com](mailto:gentech@gentechin.com); [info\\_gentech@satyam.net.in](mailto:info_gentech@satyam.net.in); Web : [www.gentechin.com](http://www.gentechin.com)

## Single Color Reagents

## Antibodies to Human Leukocytes

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD1a	MCD1a	Mo IgG1	Thymocytes	FITC	IQP-126F	1.0 mL for 100 test	CE	49,390.0
CD1a	MCD1a	Mo IgG1	Thymocytes	R-PE	IQP-126R	1.0 mL for 100 test	CE	59,500.0
CD1a	MCD1a	Mo IgG1	Thymocytes	PURE	IQP-126P	1.0 mL for 100 test		30,800.0
CD2	B-E2	Mo IgG2b	T cells, E rosette rec, thymocytes, NK cells	FITC	IQP-100F	1.0 mL for 100 test	CE	49,390.0
CD2	B-E2	Mo IgG2b	T cells, E rosette rec, thymocytes, NK cells	R-PE	IQP-100R	1.0 mL for 100 test	CE	59,500.0
CD2	B-E2	Mo IgG2b	T cells, E rosette rec, thymocytes, NK cells	PURE	IQP-100P	1.0 mL for 100 test		30,800.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	APC	IQP-519A	1.0 mL for 100 test	CE	59,500.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	CyQ	IQP-519C	1.0 mL for 100 test	CE	61,900.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	FITC	IQP-519F	1.0 mL for 100 test	CE	49,400.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	PerCP	IQP-519PC	1.0 mL for 100 test		61,900.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	PURE	IQP-519P	1.0 mL for 100 test		30,900.0
CD3	UCHT1	Mo IgG1	Thymocytes, mature T cells	R-PE	IQP-519R	1.0 mL for 100 test	CE	59,500.0
CD3	MEM-57	Mo IgG2a	Thymocytes, mature T cells	APC	IQP-173A	1.0 mL for 100 test		59,500.0
CD3	MEM-57	Mo IgG2a	Thymocytes, mature T cells	FITC	IQP-173F	1.0 mL for 100 test		49,400.0
CD3	MEM-57	Mo IgG2a	Thymocytes, mature T cells	PURE	IQP-173P	1.0 mL for 100 test		30,900.0
CD3	MEM-57	Mo IgG2a	Thymocytes, mature T cells	QCD	IQP-173D	1.0 mL for 100 test		61,900.0
CD3	MEM-57	Mo IgG2a	Thymocytes, mature T cells	R-PE	IQP-173R	1.0 mL for 100 test		59,500.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	APC	IQP-535A	1.0 mL for 100 test	CE	59,500.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	CyQ	IQP-535C	1.0 mL for 100 test	CE	61,900.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	FITC	IQP-535F	1.0 mL for 100 test	CE	49,400.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	PURE	IQP-535P	1.0 mL for 100 test		30,900.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	QCD	IQP-535D	1.0 mL for 100 test		61,900.0
CD4	Edu-2	Mo IgG2a	T helper/inducer cells, monocytes, macrophages	R-PE	IQP-535R	1.0 mL for 100 test	CE	59,500.0
CD5	MCD5	Mo IgG2b	Thymocytes, T cells, B cell subset	FITC	IQP-103F	1.0 mL for 100 test	CE	49,400.0
CD5	MCD5	Mo IgG2b	Thymocytes, T cells, B cell subset	PURE	IQP-103P	1.0 mL for 100 test		30,900.0
CD5	MCD5	Mo IgG2b	Thymocytes, T cells, B cell subset	R-PE	IQP-103R	1.0 mL for 100 test	CE	59,500.0
CD7	B-B7	Mo IgG1	T cells, NK cells	FITC	IQP-127F	1.0 mL for 100 test	CE	49,400.0
CD7	B-B7	Mo IgG1	T cells, NK cells	PURE	IQP-127P	1.0 mL for 100 test		30,900.0
CD7	B-B7	Mo IgG1	T cells, NK cells	R-PE	IQP-127R	1.0 mL for 100 test	CE	59,500.0
CD8	MCD8	Mo IgG1	T cytotoxic/suppressor cells, NK cells	APC	IQP-104A	1.0 mL for 100 test	CE	59,500.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD8	MCD8	Mo IgG1	T cytotoxic/suppressor cells, NK cells	CyQ	IQP-104C	1.0 mL for 100 test	CE	61,900.0
CD8	MCD8	Mo IgG1	T cytotoxic/suppressor cells, NK cells	FITC	IQP-104F	1.0 mL for 100 test	CE	49,400.0
CD8	MCD8	Mo IgG1	T cytotoxic/suppressor cells, NK cells	PURE	IQP-104P	1.0 mL for 100 test		30,900.0
CD8	MCD8	Mo IgG1	T cytotoxic/suppressor cells, NK cells	R-PE	IQP-104R	1.0 mL for 100 test	CE	59,500.0
CD8	MEM-31	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	APC	IQP-180A	1.0 mL for 100 test		59,500.0
CD8	MEM-31	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	FITC	IQP-180F	1.0 mL for 100 test		49,400.0
CD8	MEM-31	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	PerCP	IQP-180PC	1.0 mL for 100 test		61,900.0
CD8	MEM-32	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	PURE	IQP-180P	1.0 mL for 100 test		30,900.0
CD8	MEM-33	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	QCD	IQP-180D	1.0 mL for 100 test		61,900.0
CD8	MEM-34	Mo IgG2a	T cytotoxic/suppressor cells, NK cells	R-PE	IQP-180R	1.0 mL for 100 test		59,500.0
CD10	B-E3	Mo IgG2a	Common acute lymphoblastic leukemia antigen	FITC	IQP-105F	1.0 mL for 100 test	CE	49,400.0
CD10	B-E3	Mo IgG2a	Common acute lymphoblastic leukemia antigen	PURE	IQP-105P	1.0 mL for 100 test		30,900.0
CD10	B-E3	Mo IgG2a	Common acute lymphoblastic leukemia antigen	R-PE	IQP-105R	1.0 mL for 100 test	CE	59,500.0
CD11a	DF1524	Mo IgG2b	Lymphocytes, granulocytes, monocytes and macrophages	PURE	IQP-151P	1.0 mL for 100 test		30,900.0
CD11a	DF1524	Mo IgG2b	Lymphocytes, granulocytes, monocytes and macrophages	R-PE	IQP-151R	1.0 mL for 100 test		59,500.0
CD11b	44	Mo IgG1	T cell subset, monocytes, granulocytes, NK cells	FITC	IQP-138F	1.0 mL for 100 test	CE	49,400.0
CD11b	44	Mo IgG1	T cell subset, monocytes, granulocytes, NK cells	PURE	IQP-138P	1.0 mL for 100 test		30,900.0
CD11b	44	Mo IgG1	T cell subset, monocytes, granulocytes, NK cells	R-PE	IQP-138R	1.0 mL for 100 test	CE	59,500.0
CD11c	B-ly6	Mo IgG1	Granulocytes, monocytes, NK cells, B and T cell subset	PURE	IQP-119P	1.0 mL for 100 test		30,900.0
CD11c	B-ly6	Mo IgG1	Granulocytes, monocytes, NK cells, B and T cell subset	R-PE	IQP-119R	1.0 mL for 100 test	CE	59,500.0
CD13	22A5	Mo IgG2a	Monocytes, granulocytes, some macrophages	FITC	IQP-112F	1.0 mL for 100 test	CE	49,400.0
CD13	22A5	Mo IgG2a	Monocytes, granulocytes, some macrophages	PURE	IQP-112P	1.0 mL for 100 test		30,900.0
CD13	22A5	Mo IgG2a	Monocytes, granulocytes, some macrophages	R-PE	IQP-112R	1.0 mL for 100 test	CE	59,500.0
CD13	WM15	Mo IgG1	Monocytes, granulocytes, some macrophages	APC	IQP-531A	1.0 mL for 100 test		59,500.0
CD13	WM15	Mo IgG1	Monocytes, granulocytes, some macrophages	FITC	IQP-531F	1.0 mL for 100 test		49,400.0
CD13	WM15	Mo IgG1	Monocytes, granulocytes, some macrophages	PURE	IQP-531P	1.0 mL for 100 test		30,900.0
CD13	WM15	Mo IgG1	Monocytes, granulocytes, some macrophages	QCD	IQP-531D	1.0 mL for 100 test		61,900.0
CD13	WM15	Mo IgG1	Monocytes, granulocytes, some macrophages	R-PE	IQP-531R	1.0 mL for 100 test		59,500.0
CD14	UCHM1	Mo IgG2a	Monocytes, macrophages, granulocytes	APC	IQP-143A	1.0 mL for 100 test	CE	59,500.0
CD14	UCHM1	Mo IgG2a	Monocytes, macrophages, granulocytes	FITC	IQP-143F	1.0 mL for 100 test	CE	49,400.0
CD14	UCHM1	Mo IgG2a	Monocytes, macrophages, granulocytes	PURE	IQP-143P	1.0 mL for 100 test		30,900.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD14	UCHM1	Mo IgG2a	Monocytes, macrophages, granulocytes	R-PE	IQP-143R	1.0 mL for 100 test	CE	59,500.0
CD15	BRA-4F1	Mo IgM	Granulocytes, monocytes	FITC	IQP-129F	1.0 mL for 100 test		49,400.0
CD15	BRA-4F1	Mo IgM	Granulocytes, monocytes	PURE	IQP-129P	1.0 mL for 100 test		30,900.0
CD15	BRA-4F1	Mo IgM	Granulocytes, monocytes	R-PE	IQP-129R	1.0 mL for 100 test		59,500.0
CD16	B-E16	Mo IgG2a	NK cells, neutrophils, granulocytes, macrophages	FITC	IQP-130F	1.0 mL for 100 test		49,400.0
CD16	B-E16	Mo IgG2a	NK cells, neutrophils, granulocytes, macrophages	PURE	IQP-130P	1.0 mL for 100 test		30,900.0
CD16	B-E17	Mo IgG2a	NK cells, neutrophils, granulocytes, macrophages	R-PE	IQP-130R	1.0 mL for 100 test		59,500.0
CD19	HD37	Mo IgG1	Pan B cells	APC	IQP-515A	1.0 mL for 100 test	CE	59,500.0
CD19	HD37	Mo IgG1	Pan B cells	CyQ	IQP-515C	1.0 mL for 100 test	CE	61,900.0
CD19	HD37	Mo IgG1	Pan B cells	FITC	IQP-515F	1.0 mL for 100 test	CE	49,400.0
CD19	HD37	Mo IgG1	Pan B cells	PURE	IQP-515P	1.0 mL for 100 test		30,900.0
CD19	HD37	Mo IgG1	Pan B cells	R-PE	IQP-515R	1.0 mL for 100 test	CE	59,500.0
CD19	LT19	Mo IgG1	Pan B cells	APC	IQP-540A	1.0 mL for 100 test		59,500.0
CD19	LT19	Mo IgG1	Pan B cells	FITC	IQP-540F	1.0 mL for 100 test		49,400.0
CD19	LT19	Mo IgG1	Pan B cells	PURE	IQP-540P	1.0 mL for 100 test		30,900.0
CD19	LT19	Mo IgG1	Pan B cells	QCD	IQP-540D	1.0 mL for 100 test		61,900.0
CD19	LT19	Mo IgG1	Pan B cells	R-PE	IQP-540R	1.0 mL for 100 test		59,500.0
CD20	B-ly1	Mo IgG1	Pan B cells	APC	IQP-108A	1.0 mL for 100 test	CE	59,500.0
CD20	B-ly1	Mo IgG1	Pan B cells	FITC	IQP-108F	1.0 mL for 100 test	CE	49,400.0
CD20	B-ly1	Mo IgG1	Pan B cells	PURE	IQP-108P	1.0 mL for 100 test		30,900.0
CD20	B-ly1	Mo IgG1	Pan B cells	R-PE	IQP-108R	1.0 mL for 100 test	CE	59,500.0
CD20	LT20	Mo IgG2a	Pan B cells	APC	IQP-546A	1.0 mL for 100 test		59,500.0
CD20	LT20	Mo IgG2a	Pan B cells	FITC	IQP-546F	1.0 mL for 100 test		49,400.0
CD20	LT20	Mo IgG2a	Pan B cells	PerCP	IQP-546PC	1.0 mL for 100 test		61,900.0
CD20	LT20	Mo IgG2a	Pan B cells	PURE	IQP-546P	1.0 mL for 100 test		30,900.0
CD20	LT20	Mo IgG2a	Pan B cells	QCD	IQP-546D	1.0 mL for 100 test		61,900.0
CD20	LT20	Mo IgG2a	Pan B cells	R-PE	IQP-546R	1.0 mL for 100 test		59,500.0
CD21	B-ly4	Mo IgG1	Mature B cells, follicular dendritic cells	FITC	IQP-109F	1.0 mL for 100 test		49,400.0
CD21	B-ly4	Mo IgG1	Mature B cells, follicular dendritic cells	PURE	IQP-109P	1.0 mL for 100 test		30,900.0
CD21	B-ly4	Mo IgG1	Mature B cells, follicular dendritic cells	R-PE	IQP-109R	1.0 mL for 100 test		59,500.0
CD22	B-ly8	Mo IgG1	Pan B cells (surface and cytoplasmatic staining)	FITC	IQP-110F	1.0 mL for 100 test		49,400.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD22	B-ly8	Mo IgG1	Pan B cells (surface and cytoplasmatic staining)	PURE	IQP-110P	1.0 mL for 100 test		30,900.0
CD22	B-ly8	Mo IgG1	Pan B cells (surface and cytoplasmatic staining)	R-PE	IQP-110R	1.0 mL for 100 test		59,500.0
CD23	B-G6	Mo IgG1	B cell differentiation antigen	PURE	IQP-154P	1.0 mL for 100 test		30,900.0
CD23	B-G6	Mo IgG1	B cell differentiation antigen	R-PE	IQP-154R	1.0 mL for 100 test		59,500.0
CD25	B-B10	Mo IgG1	Activated T and B cells and activated macrophages	FITC	IQP-125F	1.0 mL for 100 test		49,400.0
CD25	B-B10	Mo IgG1	Activated T and B cells and activated macrophages	PURE	IQP-125P	1.0 mL for 100 test		30,900.0
CD25	B-B10	Mo IgG1	Activated T and B cells and activated macrophages	R-PE	IQP-125R	1.0 mL for 100 test	CE	59,500.0
CD31	MEM-05	Mo IgG1	PECAM on monocytes, platelets, granulocytes	R-PE	IQP-552R	1.0 mL for 100 test		59,500.0
CD31	NaM28-2D1	Mo IgG1	PECAM on monocytes, platelets, granulocytes	PURE	IQP-186P	1.0 mL for 100 test		30,900.0
CD31	NaM28-2D1	Mo IgG1	PECAM on monocytes, platelets, granulocytes	R-PE	IQP-186R	1.0 mL for 100 test		59,500.0
CD33	WM53	Mo IgG1	Monocytes, granulocytes	PURE	IQP-146P	1.0 mL for 100 test		30,900.0
CD33	WM53	Mo IgG1	Monocytes, granulocytes	R-PE	IQP-146R	1.0 mL for 100 test		59,500.0
CD34	581	Mo IgG1	Bone marrow progenitor cells	APC	IQP-144A	1.0 mL for 100 test	CE	59,500.0
CD34	581	Mo IgG1	Bone marrow progenitor cells	FITC	IQP-144F	1.0 mL for 100 test	CE	49,400.0
CD34	581	Mo IgG1	Bone marrow progenitor cells	PURE	IQP-144P	1.0 mL for 100 test		30,900.0
CD34	581	Mo IgG1	Bone marrow progenitor cells	R-PE	IQP-144R	1.0 mL for 100 test	CE	59,500.0
CD36	SMØ	Mo IgM	Monocytes, macrophages, platelets, weakly B cells	FITC	IQP-113F	1.0 mL for 100 test		49,400.0
CD36	SMØ	Mo IgM	Monocytes, macrophages, platelets, weakly B cells	PURE	IQP-113P	1.0 mL for 100 test		30,900.0
CD38	T16	Mo IgG1	Plasma cells, pre-B cells, immature and activated T cells	APC	IQP-132A	1.0 mL for 100 test		59,500.0
CD38	T16	Mo IgG1	Plasma cells, pre-B cells, immature and activated T cells	PURE	IQP-132P	1.0 mL for 100 test		30,900.0
CD38	T16	Mo IgG1	Plasma cells, pre-B cells, immature and activated T cells	R-PE	IQP-132R	1.0 mL for 100 test		59,500.0
CD41	NaM12-6E1	Mo IgG1	Platelets and megakaryocytes	PURE	IQP-187P	1.0 mL for 100 test		30,900.0
CD41	NaM12-6E1	Mo IgG1	Platelets and megakaryocytes	FITC	IQP-187F	1.0 mL for 100 test		49,400.0
CD41	MEM-06	Mo IgG1	Platelets and megakaryocytes	QCD	IQP-551D	1.0 mL for 100 test		61,900.0
CD41	MEM-06	Mo IgG1	Platelets and megakaryocytes	FITC	IQP-551F	1.0 mL for 100 test		49,400.0
CD43	MT1	Mo IgG1	T cells, B cells, NK cells, monocytes	FITC	IQP-133F	1.0 mL for 100 test		49,400.0
CD43	MT1	Mo IgG1	T cells, B cells, NK cells, monocytes	PURE	IQP-133P	1.0 mL for 100 test		30,900.0
CD43	MT1	Mo IgG1	T cells, B cells, NK cells, monocytes	R-PE	IQP-133R	1.0 mL for 100 test		59,500.0
CD44	MEM-85	Mo IgG2b	Epithelial, endothelial, mesothelial and mesenchymal cells	FITC	IQP-118F	1.0 mL for 100 test		49,400.0
CD44	MEM-85	Mo IgG2b	Epithelial, endothelial, mesothelial and mesenchymal cells	PURE	IQP-118P	1.0 mL for 100 test		30,900.0
CD44	MEM-85	Mo IgG2b	Epithelial, endothelial, mesothelial and mesenchymal cells	R-PE	IQP-118R	1.0 mL for 100 test		59,500.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD45	ML2	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	APC	IQP-124A	1.0 mL for 100 test	CE	59,500.0
CD45	ML2	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	CyQ	IQP-124C	1.0 mL for 100 test	CE	61,900.0
CD45	ML2	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	FITC	IQP-124F	1.0 mL for 100 test	CE	49,400.0
CD45	ML2	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	PURE	IQP-124P	1.0 mL for 100 test		30,900.0
CD45	ML2	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	R-PE	IQP-124R	1.0 mL for 100 test	CE	59,500.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	APC	IQP-549A	1.0 mL for 100 test		59,500.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	FITC	IQP-549F	1.0 mL for 100 test		49,400.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	PerCP	IQP-549PC	1.0 mL for 100 test		61,900.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	PURE	IQP-549P	1.0 mL for 100 test		30,900.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	QCD	IQP-549D	1.0 mL for 100 test		61,900.0
CD45	MEM-28	Mo IgG1	Pan leukocytes (RA, RB, RO isoforms)	R-PE	IQP-549R	1.0 mL for 100 test		59,500.0
CD45RA	MB1	Mo IgG1	B cells, T cell subset, monocytes, granulocytes	FITC	IQP-123F	1.0 mL for 100 test	CE	49,400.0
CD45RA	MB1	Mo IgG1	B cells, T cell subset, monocytes, granulocytes	PURE	IQP-123P	1.0 mL for 100 test		30,900.0
CD45RA	MB1	Mo IgG1	B cells, T cell subset, monocytes, granulocytes	R-PE	IQP-123R	1.0 mL for 100 test	CE	59,500.0
CD45RB	MT4	Mo IgG1	B cells, T cell subset, monocytes, granulocytes	PURE	IQP-140P	1.0 mL for 100 test		30,900.0
CD45RB	MT4	Mo IgG1	B cells, T cell subset, monocytes, granulocytes	R-PE	IQP-140R	1.0 mL for 100 test		59,500.0
CD45RC	MT2	Mo IgG1	B lymphocytes, T suppressor/cytotoxic cells, NK cells	FITC	IQP-117F	1.0 mL for 100 test		49,400.0
CD45RC	MT2	Mo IgG1	B lymphocytes, T suppressor/cytotoxic cells, NK cells	PURE	IQP-117P	1.0 mL for 100 test		30,900.0
CD45RO	UCHL1	Mo IgG2a	T cells, thymocytes, monocytes, macrophages	FITC	IQP-141F	1.0 mL for 100 test		49,400.0
CD45RO	UCHL1	Mo IgG2a	T cells, thymocytes, monocytes, macrophages	PURE	IQP-141P	1.0 mL for 100 test		30,900.0
CD45RO	UCHL1	Mo IgG2a	T cells, thymocytes, monocytes, macrophages	R-PE	IQP-141R	1.0 mL for 100 test		59,500.0
CD54	B-H17	Mo IgG1	ICAM-1 on lymphoid endothelium, activated lymphocytes	FITC	IQP-181F	1.0 mL for 100 test		49,400.0
CD54	B-H17	Mo IgG1	ICAM-1 on lymphoid endothelium, activated lymphocytes	PURE	IQP-181P	1.0 mL for 100 test		30,900.0
CD55	NaM16-4D3	Mo IgG1	Decay-accelerating factor, DAF	R-PE	IQP-520R	1.0 mL for 100 test		63,700.0
CD56	MOC-1	Mo IgG1	NK cells	PURE	IQP-114P	1.0 mL for 100 test		30,900.0
CD56	MOC-1	Mo IgG1	NK cells	R-PE	IQP-114R	1.0 mL for 100 test		59,500.0
CD59	NaM172-2B5	Mo IgG1	Membrane inhibitor of reactive lysis, MIRL	FITC	IQP-521F	1.0 mL for 100 test		49,400.0
CD59	NaM172-2B5	Mo IgG1	Membrane inhibitor of reactive lysis, MIRL	R-PE	IQP-521R	1.0 mL for 100 test		59,500.0
CD61	NaM28-7D6	Mo IgG1	Platelets, monocytes, macrophages, endothelial cells	FITC	IQP-188F	1.0 mL for 100 test		49,400.0
CD61	NaM28-7D6	Mo IgG1	Platelets, monocytes, macrophages, endothelial cells	PURE	IQP-188P	1.0 mL for 100 test		30,900.0
CD61	NaM28-7D6	Mo IgG1	Platelets, monocytes, macrophages, endothelial cells	R-PE	IQP-188R	1.0 mL for 100 test		59,500.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
CD64(1)	22	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	CyQ	CD64-22C	1.0 mL		61,900.0
CD64(1)	22	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	FITC	CD64-22F	1.0 mL		47,900.0
CD64(1)	22	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	PURE	CD64-22U	1.0 mL		43,400.0
CD64(1)	22	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	R-PE	CD64-22P	1.0 mL		59,000.0
CD64(1)	32.2	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	FITC	CD64-322F	1.0 mL		47,900.0
CD64(1)	32.2	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils	PURE	CD64-322U	1.0 mL		43,400.0
CD64(1)	32.2	Mo IgG1	Fc gamma receptor 1 on macrophages, monocytes and neutrophils E	R-PE	CD64-322P	1.0 mL		59,000.0
CD69	FN50	Mo IgG1	Activated lymphocytes, monocytes and platelets	FITC	IQP-553F	1.0 mL for 100 test		52,900.0
CD71	DF1513	Mo IgG1	Activated leukocytes	FITC	IQP-152F	1.0 mL for 100 test		52,900.0
CD71	DF1513	Mo IgG1	Activated leukocytes	PURE	IQP-152P	1.0 mL for 100 test		30,900.0
CD103	B-ly7	Mo IgG1	Hairy Cell Leukemia, T cells subset	FITC	IQP-111F	1.0 mL for 100 test	CE	49,400.0
CD103	B-ly7	Mo IgG1	Hairy Cell Leukemia, T cells subset	PURE	IQP-111P	1.0 mL for 100 test		30,900.0
CD103	B-ly7	Mo IgG1	Hairy Cell Leukemia, T cells subset	R-PE	IQP-111R	1.0 mL for 100 test		59,500.0
CD106	B-K9	Mo IgG1	VCAM-1. Vascular endothelium, dendritic cells, macrophages	FITC	IQP-184F	1.0 mL for 100 test		49,400.0
CD106	B-K9	Mo IgG1	VCAM-1. Vascular endothelium, dendritic cells, macrophages	PURE	IQP-184P	1.0 mL for 100 test		30,900.0
CD138	B-A38	Mo IgG1	Plasma cells	APC	IQP-153A	1.0 mL for 100 test	CE	59,500.0
CD138	B-A38	Mo IgG1	Plasma cells	CyQ	IQP-153C	1.0 mL for 100 test	CE	61,900.0
CD138	B-A38	Mo IgG1	Plasma cells	FITC	IQP-153F	1.0 mL for 100 test	CE	49,400.0
CD138	B-A38	Mo IgG1	Plasma cells	PURE	IQP-153P	1.0 mL for 100 test		30,900.0
CD138	B-A38	Mo IgG1	Plasma cells	R-PE	IQP-153R	1.0 mL for 100 test	CE	59,500.0
CD163(1)	MAC2-158	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	Biotine	CD163-158B	1.0 mL		47,900.0
CD163(1)	MAC2-158	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	FITC	CD163-158F	1.0 mL		47,900.0
CD163(1)	MAC2-158	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	PURE	CD163-158U	1.0 mL		43,400.0
CD163(1)	MAC2-158	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	R-PE	CD163-158P	1.0 mL		57,420.0
CD163(1)	MAC2-48	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	PURE	CD163-48U	1.0 mL		43,400.0
CD163(1)	MAC2-48	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	R-PE	CD163-48P	1.0 mL		59,000.0
CD163(1)	R20	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	BIOTINE	CD163-R20B	1.0 mL		47,900.0
CD163(1)	R20	Mo IgG1	Mononuclear phagocytes, monocytes and macrophages	PURE	CD163-R20U	1.0 mL		43,400.0
CD235a	NaM10-6G4	Mo IgG2a	GpA, erythroid cells and progenitors	FITC	IQP-145F	1.0 mL		49,400.0
CD235a	NaM10-6G4	Mo IgG2a	GpA, erythroid cells and progenitors	PURE	IQP-145P	1.0 mL for 100 test		30,900.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
<b>Human MHC Antigen detection</b>								
IQ HLA-B27	ABC-m3 +B7	Mo IgG2a	HLA-B27 detection with cocktail of A-B27 FITC and A-B7	PURE	FITC IQP-516F	1.0 mL for 100 test		90,300.0
HLA-B27 (2)	FD705	Mo IgG2b	HLA-B27 detection 0.5mL- 50 tests	FITC	IQP-B27F50X	0,5 mL for 50 tests	CE	67,680.0
anti-HLA-DR	BRA30	Mo IgG2a	B cells, monocytes, macrophages, activated T cells	CyQ	IQP-134C	1.0 mL for 100 test	CE	61,900.0
anti-HLA-DR	BRA30	Mo IgG2a	B cells, monocytes, macrophages, activated T cells	FITC	IQP-134F	1.0 mL for 100 test	CE	49,400.0
anti-HLA-DR	BRA30	Mo IgG2a	B cells, monocytes, macrophages, activated T cells	PURE	IQP-134P	1.0 mL for 100 test		30,900.0
anti-HLA-DR	BRA30	Mo IgG2a	B cells, monocytes, macrophages, activated T cells	R-PE	IQP-134R	1.0 mL for 100 test	CE	59,500.0
anti-HLA-DR	MEM-12	Mo IgG1	B cells, monocytes, macrophages, activated T cells	APC	IQP-550A	1.0 mL for 100 test		59,500.0
anti-HLA-DR	MEM-12	Mo IgG1	B cells, monocytes, macrophages, activated T cells	FITC	IQP-550F	1.0 mL for 100 test		49,400.0
anti-HLA-DR	MEM-12	Mo IgG1	B cells, monocytes, macrophages, activated T cells	PURE	IQP-550P	1.0 mL for 100 test		30,900.0
anti-HLA-DR	MEM-12	Mo IgG1	B cells, monocytes, macrophages, activated T cells	QCD	IQP-550D	1.0 mL for 100 test		61,900.0
anti-HLA-DR	MEM-12	Mo IgG1	B cells, monocytes, macrophages, activated T cells	R-PE	IQP-550R	1.0 mL for 100 test		59,500.0
<b>Detection of Leukocyte Enzyme-associated Molecules</b>								
anti-MPO	266.6K2	Mo IgG1	Myeloid cells, intracellular MPO	FITC	IQP-155F	1.0 mL for 100 test	CE	49,400.0
anti-MPO	266.6K2	Mo IgG1	Myeloid cells, intracellular MPO	PURE	IQP-155P	1.0 mL for 100 test		30,900.0
anti-TdT	HT-6	Mo IgG1	Lymphoid progenitor cells and precursor cells	FITC	IQP-149F	1.0 mL for 100 test		29,000.0
anti-TdT	HT-6	Mo IgG1	Lymphoid progenitor cells and precursor cells	FITC	IQP-150F	0,5 mL for 50 tests		114,200.0
<b>Markers for Cell Proliferation</b>								
anti-BrdU	Mobu-1	Mo IgG1	Specific for BdU incorporated into nuclear DNA	FITC	IQP-170F	1.0 mL for 100 test		49,400.0
anti-BrdU	Mobu-1	Mo IgG1	Specific for BdU incorporated into nuclear DNA	PURE	IQP-170P	1.0 mL for 100 test		30,900.0
anti-p53	BP53.12	Mo IgG2a	Human carcinoma	FITC	IQP-171F	1.0 mL for 100 test		49,400.0
anti-p53	BP53.12	Mo IgG2a	Human carcinoma	PURE	IQP-171P	1.0 mL for 100 test		30,900.0
<b>Human Immunoglobulin detection</b>								
anti-IgA	Polyclonal Go	Go F(ab') <sub>2</sub>	Goat anti-human IgA F(ab') <sub>2</sub>	FITC	IQP-508F	1.0 mL for 100 test		49,400.0
anti-IgD	Polyclonal Go	Go F(ab') <sub>2</sub>	Goat anti-human IgD F(ab') <sub>2</sub>	FITC	IQP-509F	1.0 mL for 100 test		49,400.0
anti-IgG	Polyclonal Go	Go F(ab') <sub>2</sub>	Goat anti-human IgG F(ab') <sub>2</sub>	FITC	IQP-510F	1.0 mL for 100 test		49,400.0
anti-IgM	Polyclonal Go	Go F(ab') <sub>2</sub>	Goat anti-human IgM F(ab') <sub>2</sub>	FITC	IQP-511F	1.0 mL for 100 test		49,400.0
<b>Diagnosis of Clonality</b>								
anti-kappa	NaM76-5F3	Mo IgG2a	50 - 70% of normal B lymphocytes	FITC	IQP-517F	1.0 mL for 100 test		49,400.0
anti-kappa	NaM76-5F3	Mo IgG2a	50 - 70% of normal B lymphocytes	R-PE	IQP-517R	1.0 mL for 100 test		59,500.0
anti-lambda	NaM79-8E6	Mo IgG1	30 - 50% of normal B lymphocytes	FITC	IQP-518F	1.0 mL for 100 test		49,400.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
anti-lambda	NaM79-8E6	Mo IgG1	30 - 50% of normal B lymphocytes	R-PE	IQP-518R			59,500.0
<b>Isotype Controls</b>								
IgG1	MCG1	Mo IgG1		APC	IQP-191A	1.0 mL for 100 test	CE	45,100.0
IgG1	MCG1	Mo IgG1		CyQ	IQP-191C	1.0 mL for 100 test	CE	45,000.0
IgG1	MCG1	Mo IgG1		FITC	IQP-191F	1.0 mL for 100 test	CE	36,300.0
IgG1	MCG1	Mo IgG1		PURE	IQP-191P	1.0 mL for 100 test		22,300.0
IgG1	MCG1	Mo IgG1		R-PE	IQP-191R	1.0 mL for 100 test	CE	42,300.0
IgG2a	MCG2a	Mo IgG2a		APC	IQP-192A	1.0 mL for 100 test	CE	42,300.0
IgG2a	MCG2a	Mo IgG2a		CyQ	IQP-192C	1.0 mL for 100 test		45,000.0
IgG2a	MCG2a	Mo IgG2a		FITC	IQP-192F	1.0 mL for 100 test	CE	36,300.0
IgG2a	MCG2a	Mo IgG2a		PURE	IQP-192P	1.0 mL for 100 test		22,300.0
IgG2a	MCG2a	Mo IgG2a		R-PE	IQP-192R	1.0 mL for 100 test	CE	42,300.0
IgG2b	MCG2b	Mo IgG2b		FITC	IQP-193F	1.0 mL for 100 test	CE	36,300.0
IgG2b	MCG2b	Mo IgG2b		PURE	IQP-193P	1.0 mL for 100 test		22,300.0
IgG2b	MCG2b	Mo IgG2b		R-PE	IQP-193R	1.0 mL for 100 test	CE	42,300.0
IgM	MCGM	Mo IgM		FITC	IQP-194F	1.0 mL for 100 test	CE	36,300.0
IgM	MCGM	Mo IgM		PURE	IQP-194P	1.0 mL for 100 test		22,300.0
IgM	MCGM	Mo IgM		R-PE	IQP-194R	1.0 mL for 100 test	CE	42,300.0
Ig	Polyclonal Rb	Rb Ig		FITC	IQP-195F	1.0 mL for 100 test		36,300.0
Ig	Polyclonal Go	Go Ig		FITC	IQP-196F	1.0 mL for 100 test		36,300.0
Ig	Polyclonal Go	Go Ig		R-PE	IQP-196R	1.0 mL for 100 test		42,300.0
2nd Step Reagent	Polyclonal Rb	Rb F(ab') <sub>2</sub>	Rabbit anti-mouse IgG F(ab') <sub>2</sub>	FITC	IQP-190F	0.5 mL (10x) - 100 tests		22,000.0
2nd Step Reagent	Polyclonal Rb	Rb F(ab') <sub>2</sub>	Rabbit anti-mouse IgG F(ab') <sub>2</sub>	R-PE	IQP-190R	0.5 mL (10x) - 100 tests		30,000.0
<b>Single Color Reagents Against Human Cytokines</b>								
IFN-g	45-14	Mo IgG1	T cells, NK cells	APC	IQP-160A	0.5 mL for 50 tests		39,400.0
IFN-g	45-14	Mo IgG1	T cells, NK cells	FITC	IQP-160F	0.5 mL for 50 tests		33,000.0
IFN-g	45-14	Mo IgG1	T cells, NK cells	PURE	IQP-160P	0.5 mL for 50 tests		23,700.0
IFN-g	45-14	Mo IgG1	T cells, NK cells	R-PE	IQP-160R	0.5 mL for 50 tests		38,300.0
IL-1 $\beta$	B-A15	Mo IgG1	Monocytes, macrophages	PURE	IQP-167P	0.5 mL for 50 tests		24,000.0
IL-1 $\beta$	B-A15	Mo IgG1	Monocytes, macrophages	R-PE	IQP-167R	0.5 mL for 50 tests		38,300.0
IL-2	N7-48A	Mo IgG2a	T helper 1 cells	PURE	IQP-161P	0.5 mL for 50 tests		24,000.0

Item	clone	Isotype	Cellular expression	label	product code	unit size	CE	Price Rs
IL-2	N7-48A	Mo IgG2a	T helper 1 cells	R-PE	IQP-161R	0.5 mL for 50 tests		38,300.0
IL-4	8F-12	Mo IgG1	T cells, mast cells and stromal cells	PURE	IQP-162P	0.5 mL for 50 tests		24,000.0
IL-4	8F-12	Mo IgG1	T cells, mast cells and stromal cells	R-PE	IQP-162R	0.5 mL for 50 tests		38,300.0
IL-6	B-E8	Mo IgG1	Monocytes, macrophages, fibroblasts and endothelial, T and B cells	PURE	IQP-164P	0.5 mL for 50 tests		24,000.0
IL-6	B-E8	Mo IgG1	Monocytes, macrophages, fibroblasts and endothelial, T and B cells	R-PE	IQP-164R	0.5 mL for 50 tests		38,300.0
IL-10	B-N10	Mo IgG1	T helper 2 cells	FITC	IQP-165F	0.5 mL for 50 tests		33,300.0
IL-10	B-N10	Mo IgG1	T helper 2 cells	PURE	IQP-165P	0.5 mL for 50 tests		24,000.0
IL-10	B-N10	Mo IgG1	T helper 2 cells	R-PE	IQP-165R	0.5 mL for 50 tests		38,300.0
IL-12	B-P24	Mo IgG1	Antigen presenting cells e.g. monocytes, B cells and dendritic cells	PURE	IQP-168P	0.5 mL for 50 tests		24,000.0
IL-12	B-P24	Mo IgG1	Antigen presenting cells e.g. monocytes, B cells and dendritic cells	R-PE	IQP-168R	0.5 mL for 50 tests		39,000.0
IL-13	B-B13	Mo IgG1	Activated T helper 2 cells	PURE	IQP-166P	0.5 mL for 50 tests		24,000.0
IL-13	B-B13	Mo IgG1	Activated T helper 2 cells	R-PE	IQP-166R	0.5 mL for 50 tests		38,300.0
TGF- 1	TB21	Mo IgG1	Many cells including T cells and monocytes	PURE	IQP-169P	0.5 mL for 50 tests		45,100.0
TGF- 1	TB21	Mo IgG1	Many cells including T cells and monocytes	R-PE	IQP-169R	0.5 mL for 50 tests		71,900.0
TNF-	B-C7	Mo IgG1	T cells, B cells and monocytes	PURE	IQP-163P	0.5 mL for 50 tests		24,000.0
TNF-	B-C7	Mo IgG1	T cells, B cells and monocytes	R-PE	IQP-163R	0.5 mL for 50 tests		38,300.0
<b>Single Color Reagents Against Human Erythrocyte Antigens</b>								
CD235a	NaM10-6G4	Mo IgG2a	GpA, erythroid cells and progenitors	FITC	IQP-145F	1.0 mL for 100 test		49,400.0
CD235a	NaM10-6G4	Mo IgG2a	GpA, erythroid cells and progenitors	PURE	IQP-145P	1.0 mL for 100 test		30,900.0
Anti-D	NaTH109-1G2	Hu IgG	Rh D antigen (RH1)	R-PE	IQP-513R	1.0 mL for 100 test		59,500.0

Go = Goat  
Ha = Hamster  
Hu = Human  
Mo = Mouse  
Rb = Rabbit  
Rt = Rat

- (1) Exclusively distributed by IQ Products outside the US for Trillium Diagnostics, LLC
- (2) Product available in the Netherlands only
- (3) CyQ label. Excitation: 488 nm - Emission: 670 nm. 4 CyQ is a tandem dye of PE and Cy5.18
- (4) QCD label. Excitation: 488 nm - Emission: 599 nm. 6 QCD is a tandem dye of PE and DY590 (Dyomics590), spectral similar to PE-Texas Red (ECD) and PE-Alexa Fluor®590. DY590 has spectral similar to Texas Red and Alexa Fluor®590.
- (5) RQ label. Excitation: 488 nm - Emission: 782 nm. RQ is a tandem dye of PE and DY747(Dyomics 747), spectral similar to PE-Cy®7. DY747 is a heptamethine analogue of Dyomics647 (a pentamethine spectral similar to Cy®5 and Alexa Fluor®647).
- (6) FRQ label. Excitation: 595-640 nm - Emission: 782 nm. 8 FRQ is a tandem dye of 5 APC and DY747(Dyomics 747), spectral similar to 5 APC-Cy®7. DY747 is a heptamethine analogue of Dyomics647 (a pentamethine spectral similar to Cy®5 and Alexa Fluor®647).

## Dual Color Reagents

FITC	clone	R-PE	clone	Product code	unit size	CE	Price Rs
CD2	B-E2	CD19	HD37	IQP-252FR	1 mL for 50 tests		46,900.0
CD3	UCHT1	CD4	Edu-2	IQP-253FR	1 mL for 50 tests	CE	46,900.0
CD3	UCHT1	CD8	MCD8	IQP-254FR	1 mL for 50 tests	CE	46,900.0
CD3	UCHT1	CD16	B-E16	IQP-256FR	1 mL for 50 tests		46,900.0
		CD56	MOC-1				46,900.0
CD3	UCHT1	CD19	HD37	IQP-249FR	1 mL for 50 tests	CE	46,900.0
CD3	UCHT1	CD25	B-B10	IQP-257FR	1 mL for 50 tests		46,900.0
CD3	UCHT1	HLA-DR	BRA30	IQP-258FR	1 mL for 50 tests		46,900.0
CD4	Edu-2	CD8	MCD8	IQP-259FR	1 mL for 50 tests	CE	46,900.0
CD4	Edu-2	CD25	B-B10	IQP-260FR	1 mL for 50 tests		46,900.0
CD4	Edu-2	CD45RA	MB1	IQP-261FR	1 mL for 50 tests		46,900.0
CD5	MCD5	CD19	HD37	IQP-262FR	1 mL for 50 tests		46,900.0
CD5	MCD5	CD20	B-ly1	IQP-210FR	1 mL for 50 tests		46,900.0
CD7	B-B7	CD56	MOC-1	IQP-241FR	1 mL for 50 tests		46,900.0
CD8	MCD8	CD25	B-B10	IQP-213FR	1 mL for 50 tests		46,900.0
CD8	MCD8	CD38	T16	IQP-230FR	1 mL for 50 tests		46,900.0
CD8	MCD8	CD56	MOC-1	IQP-214FR	1 mL for 50 tests		46,900.0
CD10	B-E3	CD19	HD37	IQP-264FR	1 mL for 50 tests	CE	46,900.0
CD16	B-E16	CD56	MOC-1	IQP-221FR	1 mL for 50 tests		46,900.0
CD19	HD37	CD10	B-E3	IQP-265FR	1 mL for 50 tests		46,900.0
CD22	B-ly8	CD5	MCD5	IQP-242FR	1 mL for 50 tests		46,900.0
CD22	B-ly8	CD56	MOC-1	IQP-244FR	1 mL for 50 tests		46,900.0
CD45	ML2	CD14	UCHM1	IQP-228FR	1 mL for 50 tests	CE	46,900.0
CD45RA	MB1	CD4	Edu-2	IQP-266FR	1 mL for 50 tests		46,900.0
CD45RO	UCHL1	CD4	Edu-2	IQP-267FR	1 mL for 50 tests		46,900.0
CD59	NaM172-2B5	CD55	NaM16-4D3	IQP-273FR	1 mL for 50 tests		46,900.0
CD103	B-ly7	CD19	HD37	IQP-268FR	1 mL for 50 tests	CE	46,900.0
CD103	B-ly7	CD22	B-ly8	IQP-227FR	1 mL for 50 tests	CE	46,900.0
CD138	B-A38	CD38	T16	IQP-243FR	1 mL for 50 tests		46,900.0
CD138	B-A38	CD56	MOC-1	IQP-236FR	1 mL for 50 tests		46,900.0

FITC	clone	R-PE	clone	Product code	unit size	CE	Price Rs
HLA-DR	BRA30	CD3	UCHT1	IQP-269FR	1 mL for 50 tests		46,900.00
HLA-DR	BRA30	CD8	MCD8	IQP-217FR	1 mL for 50 tests		46,900.00
anti-Kappa	NaM76-5F3	CD19	HD37	IQP-246FR	1 mL for 50 tests		46,900.00
anti-Lambda	NaM79-8E6	CD19	HD37	IQP-247FR	1 mL for 50 tests		46,900.00
anti-Kappa	NaM76-5F3	anti-Lambda	NaM79-8E6	IQP-248FR	1 mL for 50 tests		46,900.00

### Isotype Controls

IgG1	MCG1	IgG1	MCG1	IQP-291FR	1 mL for 50 tests	CE	33,000.0
IgG1	MCG1	IgG2a	MCG2a	IQP-292FR	1 mL for 50 tests	CE	33,000.0
IgG2b	MCG2b	IgG1	MCG1	IQP-293FR	1 mL for 50 tests		33,000.0
IgG2a	MCG2a	IgG1	MCG1	IQP-294FR	1 mL for 50 tests	CE	33,000.0
IgG2a	MCG2a	IgG2b	MCG2b	IQP-276FR	1 mL for 50 tests		33,000.0
IgG1	MCG1	IgG2b	MCG2b	IQP-277FR	1 mL for 50 tests		33,000.0
IgG1	MCG1	IgG2a	MCG2a	IQP-278FR	1 mL for 50 tests		33,000.0
		IgG1	MCG1				33,000.0
IgG1	MCG1	IgG2a	MCG2a	IQP-279FR	1 mL for 50 tests		33,000.0
		IgG2b	MCG2b		1 mL for 50 tests		33,000.0
Polyclonal Ig	Rabbit	Polyclonal Ig	Goat	IQP-295FR	1 mL for 50 tests		33,000.0
Polyclonal Ig	Rabbit	IgG1	MCG1	IQP-296FR	1 mL for 50 tests		33,000.0

FITC	Clone	R-PE	Clone	CyQ	Clone	Product code	Unit size	CE	Price Rs
CD3	UCHT1	CD4	Edu-2	CD45	ML2	IQP-417FRC	1 mL for 50 tests		66,000.0
CD3	UCHT1	CD8	MCD8	CD45	ML2	IQP-418FRC	1 mL for 50 tests		66,000.0
CD3	UCHT1	CD19	HD37	CD45	ML2	IQP-419FRC	1 mL for 50 tests		66,000.0
CD4	Edu-2	CD8	MCD8	CD3	UCHT1	IQP-420FRC	1 mL for 50 tests	CE	66,000.0
CD15	BRA-4F1	CD14	UCHM1	CD45	ML2	IQP-402FRC	1 mL for 50 tests		66,000.0
CD16	B-E16	CD19	HD37	CD3	UCHT1	IQP-423FRC	1 mL for 50 tests		66,000.0
anti-Kappa	NaM76-5F3	anti-Lambda	NaM79-8E6	CD19	HD37	IQP-411FRC	1 mL for 50 tests		66,000.0
anti-Kappa	NaM76-5F3	anti-Lambda	NaM79-8E6	CD138	B-A38	IQP-416FRC	1 mL for 50 tests		66,000.0
IgG1	MCG1	IgG1	MCG1	IgG1	MCG1	IQP-491FRC	1 mL for 50 tests		41,600.0
IgG1	MCG1	IgG2a	MCG2a	IgG1	MCG1	IQP-425FRC	1 mL for 50 tests		41,600.0
IgG2a	MCG2a	IgG1	MCG1	IgG1	MCG1	IQP-426FRC	1 mL for 50 tests	CE	41,600.0
Polyclonal Ig	Rabbit	Polyclonal Ig	Goat	IgG1	MCG1	IQP-492FRC	1 mL for 50 tests		41,600.0

**Triple Color Reagents****Isotype Controls****Transplant related products**

Item	Description	Product Code	Unit Size	CE	Price Rs
CMV Brite Turbo kit	Complete kit of the diagnosis of an active CMV infection (2 hours)	VIR-CMV 110	110	CE	126,000.0
CMV Brite kit	Complete kit of the diagnosis of an active CMV infection (4 hours)	VIR-CMV 40	40 tests	CE	87,300.0
CMV Brite kit	Complete kit of the diagnosis of an active CMV infection (4 hours)	VIR-CMV 100	100 tests	CE	149,000.0
CMV C10/C11	Cocktail Mix of 2 CMV pp65 monoclonal antibodies	VIR-CMV C10/C11	200 tests	CE	58,200.0
FITC Conjugate	Sheep anti-mouse-FITC	VIR-FITC	200 tests	CE	28,500.0
CMV Control slides	CMV pp65 positive and negative cytopots	VIR-CMV CS05	set of 5	CE	15,500.0

**Transfer prices of kit components****CMV Brite Turbo kit**

Reagent	Volume	Price Rs
Reagent A	200 ml	11,400.0
Reagent B	290 ml	12,450.0
Reagent C	290 ml	12,450.0
Reagent D	4 ml	11,400.0
Reagent E	4 ml	11,400.0

**CMV Brite kit 100 test**

Reagent	Volume	Price Rs
Reagent A	100 ml	9,500.0
Reagent B	30 ml	7,200.0
Reagent C	290 ml	12,450.0
Reagent D	290 ml	12,450.0
Reagent E	15 ml	11,400.0
Reagent F	5.9 ml	12,450.0
Reagent G	5.9 ml	12,450.0

**CMV Brite kit 40 test**

Reagent	Volume	Price Rs
Reagent A	50 ml	7,200.0
Reagent B	15 ml	5,800.0
Reagent C	290 ml	12,450.0
Reagent D	290 ml	12,450.0
Reagent E	15 ml	11,400.0
Reagent F	3.4 ml	9,500.0
Reagent G	3.4 ml	9,500.0

**Miscellaneous Reagents****Detection of Apoptosis**

Item	Cellular Expression	Label	Product Code	Unit Size	CE	Price Rs
Annexin V	Binds phosphatidyl serine	APC	IQP-120A	1 ml for 100 Tests		59,500.0
Annexin V	Binds phosphatidyl serine	FITC	IQP-120F	1 ml for 100 Tests		49,400.0
Annexin V	Binds phosphatidyl serine	R-PE	IQP-120R	1 ml for 100 Tests		59,500.0
Propidium iodide	Stains dead and necrotic cells	PURE	IQP-121	1 ml for 100 Tests		14,000.0
Merocyanine 540	Binds membranes of apoptotic cells	PURE	IQP-198	1 ml for 100 Tests		14,000.0
Phosphatidyl Serine Detection Kit	Complete kit containing: Annexin V, propidium iodide and calcium buffer	FITC	IQP-116F	kit for 100 test		54,000.0
Phosphatidyl Serine Detection Kit	Complete kit containing: Annexin V, propidium iodide and calcium buffer	R-PE	IQP-116R	kit for 100 test		61,900.0

Absolute counting					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Bright Count Microspheres™	Single platform system for absolute counting	IQP-391	50 tests		33,900.0
Bright Count Microspheres™	Single platform system for absolute counting	IQP-392	100 tests		58,100.0
Lysing and Permeabilization					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
IQ Lyse	'Ready to use' erythrocyte lysing solution (3)	IQP-199	250 test		16,650.0
IQ Brightlyse	'Ready to use' erythrocyte lysing solution designed for use with Bright Count Microspheres™	IQP-393	250 test		37,800.0
IQ Starfiqs	Intracellular staining fixation and permeabilization solution	IQP-200	50 test		16,650.0
Hybridokine™					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Hybridokine™	Culture Medium additive	Hykine 12.5	12.5 mL		12,100.0
Kits for the Detection of Human Cytokines					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Cydetect™ kit	Kit containing all reagents needed for stimulation and permeabilization	IQP-366	Kit for 100 Tests		72,900.0
Kits for the Detection of Human Fetal Erythrocytes					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Fetal Cell Count™ kit	Complete assay for routine diagnosis of fet-maternal transfusion (FMT) using anti HbF and anti-CA	IQP-379	Kit for 25 Tests	CE	59,700.0
Fetalrol	Tri-level stabilized blood controls with known human fetal erythrocytes content in human adult blood (1)	IQP-370 FT	2x2 vials of each level	CE	32,200.0
Kits for the Detection of Sepsis					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Leuko64™	Flow cytometry assay for the detection of inflammation and tissue injury by leukocyte CD64 measurement (1)	LK-064-75	Kit for 75 Tests	CE	88,800.0
Leuko64™	Flow cytometry assay for the detection of inflammation and tissue injury by leukocyte CD64 measurement (1)	LK-064-250	Kit for 250 Tests	CE	251,000.0
Leuko64™ - H	Leukocyte CD64 measurement assay for the Abbott Cell-Dyn 4000™ or Sapphire™ hematology analyser (1)	LK-64H-100	Kit for 100 Tests	CE	222,000.0
Neopterin ELISA Kit	Neopterin ELISA assay quantification of neopterin in serum, plasma and urine	IQP-386	Kit for 96 Tests	CE	59,500.0
Kits for Haematology					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
HIT Alert kit	A complete kit for the Reliable Diagnosis of Heparin-Induced Thrombocytopenia (HIT)	IQP-396	Kit for 25 test		45,000.0
Kits for Inflammation					
Item	Cellular Expression	Product Code	Unit Size	CE	Price Rs
Macro163™	sCD163 Elisa assay for macrophage and monocyte activation	IQP-383	Kit for 96 Tests		79,900.0
sTREM-1 ELISA	sTREM-1 Elisa assay quantification of soluble Trem-1 in serum and plasma	IQP-385	Kit for 96 Tests		79,900.0



## Ordering Information

Place your order from anywhere in India by Telephone, Fax or E-mail  
Call our customer service reps at Ph : 91-11-2644 5263, 2642 9843

---

### Mail Your Order to:

Gentech Diagnostics P. Ltd.  
234-A, First Floor, Chandi Plaza, Sant Nagar,  
East of Kailash, New Delhi - 110065 India

---

### Fax Your Order to:

91-11-2644 5269

---

**E-mail : [info\\_gentech@satyam.net.in](mailto:info_gentech@satyam.net.in); [gentech@gentechin.com](mailto:gentech@gentechin.com)**

---

**HELP** us to speed up your order by clearly stating the following information on the Order Form attached with the Price-List or on Telephone.

1. Your name and contact no.
  2. Complete Billing & Shipping address.
  3. Purchase Order Number, (no order will be shipped without official P.O. #)
  4. Catalog number, Description and Quantities.
- 

### Minimum Order

There are no minimum order requirements on orders. All prices are MRP as stated in Catalog.

---

### Delivery

We endeavor 7-14 days delivery anywhere in India, subject to availability of the products. We maintain the inventory of Products of different storage conditions and make our best efforts to deliver most of your requirements at a shortest possible time. But sometimes, due to backorders with our overseas supplier we may not be able to meet the delivery schedule.

---

### Dry-Ice Shipment

Minimum packing charges of Rs. 500/- will be charged irrespective of the volume of order.

---

### Taxes

Central Sales Tax/VAT will be charged extra. If you are approved by the Govt. of India for exemption for products used for education or research purpose, please submit the necessary Declaration Form (Form "C" or "D" etc). In the absence of these forms accompanying your order, we will be charging applicable taxes. Octroi and other levies if any will be to the buyers account wherever applicable.

---

### Warranty

The Products under the price list meet rigorous quality standards and are warranted by our manufacturing supplier. Gentech warranty all Products and other supplies to be free from any defect in material or workmanship to the purchaser. For a period of three months from the date of invoice, Gentech will correct either by repair or replacement at the company's sole discretion, any defect in material or workmanship without charge excluding any misuse, mishandling, abnormal storage conditions, accident or act of God or modification of the nature of the Product, in its sole discretion. Without prior authorization Gentech does not cover the liability. Determination of whether any product has been subject to misuse of abuse will be made solely by the company upon being convinced of the genuineness of such claims, if any made by the user.

---

### Damage Claims

The shipment should be inspected immediately on arrival. Damaged goods should be reported to the carrier or our local representative or our dealer within 24 hours and the damage report should be sent to the Gentech head office.

## General Terms & Conditions for Sale:-

The formal Purchase Orders are accepted by Gentech Diagnostics P. Ltd. subject to terms and conditions stated hereinafter.

- 1) We will deliver technically & commercially clear orders within 6-8 weeks of their receipt, unless the product (s) are on back-order or discontinued by the principal manufacturer. Products / Services may get deleted between the proforma & the confirmed purchase order which may or may not be served with alternatives, as to be confirmed through the Order Confirmation. The purchasing authority is to study the order confirmation by Gentech and intimate any variances before actual shipment to take place, per shipping schedule given on the Order Confirmation.
- 2) Part shipment are generally considered approved by Gentech, unless specifically stated on your purchase order or intimated to be unacceptable through our Order Confirmation.
- 3) The Prices shown in the Price List are ex-Delhi, Warehouse prices.
- 4) Payment terms: All supplies are subject to advance payment receipts (unless, otherwise stated in the proforma or Order Confirmation). Freight, Handling and Insurance shall apply separately, unless included already. Any supplies made on sight draft that remain outstanding for more than 30 days shall attract interest @2.5% per month from the date of invoice upto the date of payment, calculated on the week basis, taking part-week as full. The Purchaser on the pretext of part shipment so confirmed through the Order Confirmation is not to hold back the Payment or component thereof. Payment should be made for each such Part shipment in full.
- 5) Gentech does not indemnify the Purchaser for any loss, injury, or damage of any kind, direct or indirect, howsoever caused. The maximum liability of Gentech is not greater than the Product value supplied, for which the Purchaser has to provide all necessary details to satisfy Gentech to receive any claims. Any Purchaser claim(s) should be reported within Three Days from the receipt of the stocks, or else such claims shall not be entertained.
- 6) Gentech shall have lien on the Products for the unpaid amount. Gentech through its authorized representative shall be entitled to do such things as may be necessary to repossess the Products.
- 7) Any variances made in quoted products; quantities, prices, terms or conditions than those given in the proforma at the time of placing the order, shall not be a binding on Gentech, unless agreed upon in writing.
- 8) Prices are subject to change from time to time. Please ask for Proforma before placing the Confirmed Order.
- 9) The Products or Services are for laboratory research use only and are not intended for human or animal diagnostic, therapeutic, or other clinical unless otherwise stated.
- 10) The Purchaser shall pay for all legal Taxes / Octroi / Levies etc. as applicable from time to time. Documents pertaining to any Tax concession/way bills etc should accompany the order.
- 11) Certain Products are sold in a minimum quantity. Please inquire at the time of placing the order.

**Important:** Please refer to these Terms and Conditions page before ordering. The same is requested on each product Segment Page inside the Price List.

**For Rate Contract Customers:** Some of the above referred Terms & Conditions could be different in your case which would be in Rate Contract proposals. However the remaining above referred Terms and Conditions shall still apply in addition.